


# PCR Reagents





### ***TaqNova* Polymerase Family:**

- Validated for a wide range of templates and applications
- Extreme yield even with minimal amounts of enzyme and low DNA concentrations
- Available also for direct gel loading of PCR products, for GC-rich templates and as ready-to-use all-in one tube master mixes

# PCR REAGENTS

We provide a wide range of reagents used in various PCR techniques, including:

- Thermostable DNA polymerases (**TaqNova**, **Hypernova**, **LongNova**) along with reaction buffers and magnesium salts for reaction optimization,
- **dNTPs** available as a equimolar mixture of **dATP**, **dCTP**, **dGTP** and **dTTP** or as a set of four separate single-nucleotide solutions,
- various PCR enhancers (**PCR Anti-inhibitor**, **TaqSSB protein**, **5x GC-Additive** etc.),
- ready-to-use PCR master mixes containing all the necessary PCR reagents in a single tube.

We also offer design and/or synthesis of PCR primers and molecular probes for the Real-Time PCR.

	PRODUCT	5' → 3' exonuclease activity	3' → 5' exonuclease activity	Taq polymerase	Pyrococcus polymerase	Proofreading	Mixes	Diagnostic PCR	Long Range PCR (up to 20 kb)	Blunt End Cloning	TA Cloning	GC-rich templates	Difficult templates	Multiplex PCR	Site-directed mutagenesis	Direct load
POLYMERASES	<i>TaqNova</i>	✓		✓			✓	✓			✓			✓		
	<i>TaqNova-RED</i>	✓		✓			✓	✓			✓			✓		✓
	<i>TaqNovaGC</i>	✓		✓				✓			✓	✓	✓			
	<i>Hypernova</i>		✓		✓	✓	✓	✓		✓			✓	✓	✓	
	<i>Hypernova-RED</i>		✓		✓	✓	✓	✓		✓			✓	✓	✓	✓
	<i>LongNova</i>	✓	✓	✓	✓	✓	✓		✓		✓			✓	✓	
	<i>LongNova-RED</i>	✓	✓	✓	✓	✓	✓		✓		✓			✓	✓	✓
ENHANCERS	<i>PCR Anti-inhibitor</i>												✓	✓		
	<i>TaqSSB Protein</i>											✓	✓	✓		
	<i>5x GC-Additive</i>											✓	✓	✓		

## TaqNova Polymerase

*TaqNova* DNA Polymerase is a 94 kDa recombinant, thermostable *Taq* DNA polymerase isolated from *Thermus aquaticus*. It is recommended for a wide range of applications, which require DNA synthesis in extremely high temperatures. The *TaqNova* polymerase is a universal and easy-to-use DNA polymerase, which works rapidly and effectively in various PCR conditions. The enzyme catalyses DNA synthesis in a 5'→3' directions, it does not show a 3'→5' exonuclease activity, however it has a 5'→3' exonuclease activity.

- Consistent results
- Suitable for a wide range of applications
- Recombinant enzyme of high purity
- Extreme yield with minimal amounts of enzyme and little optimisation
- Half-life of the enzyme is 45 minutes at 94°C
- Amplifies fragments of up to 5 kb
- Leaves 'A' overhangs

### Applications:

- Efficient amplification of short and medium size DNA sequences
- Diagnostic PCR

## Quick and simple PCR set-up

We propose **TaqNova-RED** DNA Polymerase and **2x PCR TaqNova-RED** to simplify PCR operations. This is vital, especially for routine applications e.g. diagnostic assays.


### TaqNova-RED DNA Polymerase

The *TaqNova-RED* DNA Polymerase consists of *TaqNova* DNA polymerase with an inert red dye to facilitate accurate low volume pipetting and as an indicator of an enzyme addition. This dye has no adverse effect on the outcome of PCR; yields are the same as with standard *TaqNova* DNA polymerase. Use of the *TaqNova-RED* DNA Polymerase decreases the risk of making a mistake during a reaction set-up e.g. skipping the polymerase, inaccurate reagents mixing. Additionally a PCR product can be applied directly onto a gel after amplification without mixing with loading buffer.

- **Facilitate PCR reaction set-up**
- **Reduced risk of mistake, thanks to the inert red dye**
- **Direct gel loading**

### 2x PCR TaqNova-RED – PCR Mix

The **2x PCR TaqNova-RED** contains *TaqNova* DNA polymerase, dNTPs, optimal PCR buffer and all other components required for PCR, except for DNA template and primers. The **2x PCR TaqNova-RED** application reduces the number of pipetting steps and contamination risk, facilitates greater efficiency, throughput and reproducibility. The PCR mixtures prepared with the **2x PCR TaqNova-RED** can be applied onto a gel directly after amplification, no loading buffer use needed.



**Fig. 1 Efficiency and sensitivity of Taq polymerases.**

A 630 bp fragment of *ccr5* gene was amplified using the *TaqNova* polymerase and results were compared with the results obtained in a parallel reaction using *Taq* polymerases from supplier B, supplier F and supplier A. The process used serial dilutions of human genomic DNA, incubated for 2 min at 95°C, followed by 40 cycles of 30 sec at 95°C, 30 sec at 54°C, 30 sec at 72°C and a final extension 2 min at 72°C. The DNA markers used were M100-1000 and M600-1000. *TaqNova* delivers higher or comparable yields and is more or equally sensitive comparing to 3 competing products.

It saves time and facilitates agarose or polyacrylamide gel electrophoresis preparation.

- **Facilitate PCR reaction set-up**
- **Reduced risk of contamination**
- **Direct gel loading**

Ordering information			
PRODUCT	CONCENTRATION	VOLUME	CAT. NO.
<i>TaqNova</i> DNA Polymerase	2 U/μl	200 U	RP702
		500 U	RP705
		1000 U	RP710
		2500 U	RP725
	5 U/μl	200 U	RP702A
		500 U	RP705A
1000 U		RP710A	
<i>TaqNova-RED</i> DNA Polymerase	1 U/μl	200 U	RP20R
		1000 U	RP100R
<b>2x PCR TaqNova-RED</b>	0.04 U/μl	100 reactions	RP85T
		1000 reactions	RP85T-10
<i>TaqNovaGC</i> DNA Polymerase	5 U/μl	200 U	RP70-020
		1000 U	RP70-100

## GC-rich templates

Successful amplification depends of GC content and the complexity of the DNA template. The **TaqNovaGC** DNA Polymerase is an ideal tool for amplification of GC-rich templates. Special Buffer – **5x GC-Additive** changes DNA behaviour upon heating and can be used with primer – template pairs with GC-rich content which do not work well with standard PCR conditions.

The **5x GC-Additive** reduces the number of secondary structures and enables specific hybridization of primers.

- **GC-rich templates**
- **High specificity**
- **Ideal for problematic templates, which fail with standard Taq DNA polymerases**

## Hypernova Polymerase

*Hypernova* DNA Polymerase is a unique blend of a highly thermostable and proofreading modified DNA polymerase *Pwo* isolated from the hyperthermophilic archaeon *Pyrococcus woesei* and enzymes increasing yield and performance of the PCR reaction. The enzyme mix can generate very long amplicons (over 10 kb). Such long amplicons are often difficult to generate with a single DNA polymerase. *Hypernova* is a versatile and easy-to-use polymerase, since it works with many different protocols and requires minimal time consuming optimisation. The polymerase produces higher yields than most commercially available enzymes and is ideally suited for difficult PCR templates. The polymerase maintains 95% activity even after 40 cycles consisting of three one minute steps each.

### Advantages:

- High processivity (for very long amplicons)
- High yield with minimal amounts of enzyme and little optimisation
- High fidelity (proofreading activity)
- Mistake-proof during multiplex PCR
- Very specific and sensitive
- More thermostable than *Taq* polymerase

## Quick and simple PCR set-up

We propose **Hypernova-RED** DNA Polymerase and **2x PCR Hypernova-RED** to simplify PCR operations, especially vital for routine applications e.g. diagnostic assays.

### Hypernova-RED DNA Polymerase

*Hypernova-RED* DNA Polymerase consists of *Hypernova* DNA polymerase with an inert red dye to facilitate accurate low volume pipetting and as an indicator of an enzyme addition. This dye has no adverse effect on the outcome of PCR; yields are the same as with standard *Hypernova* DNA polymerase. Use of the *Hypernova-RED* DNA Polymerase decreases the risk of making a mistake during a reaction set-up e.g. skipping the polymerase, inaccurate reagents mixing. Additionally PCR product can be applied directly onto gel after amplification without mixing with loading buffer.

- **Facilitate PCR reaction set-up**
- **Reduced risk of mistake thanks to inert red dye**
- **Direct gel loading**

Ordering information			
PRODUCT	CONCENTRATION	VOLUME	CAT. NO.
<b>Hypernova</b> DNA Polymerase	2 U/μl	200 U	RP232
		1000 U	RP235
<b>Hypernova-RED</b> DNA Polymerase	1 U/μl	200 U	RP232R
		1000 U	RP235R
<b>2x PCR</b> <b>Hypernova-RED</b>	0.04 U/μl	100 reactions	RP85
		1000 reactions	RP85-10

### Applications:

- Long range PCR
- High fidelity for purposes of blunt-end PCR cloning, site-directed mutagenesis, etc.
- Multiplex PCR

### 2x PCR Hypernova-RED – PCR Mix

The **2x PCR Hypernova-RED** contains *Hypernova* DNA polymerase, dNTPs, optimal PCR buffer and all other components required for PCR, except for DNA template and primers. The **2x PCR Hypernova-RED** application reduces the number of pipetting steps and contamination risk, facilitates greater efficiency, throughput and reproducibility. The PCR mixtures prepared with **2x PCR Hypernova-RED** can be applied onto a gel directly after amplification, no loading buffer use needed. It saves time and facilitates the agarose or polyacrylamide gel electrophoresis preparation.

- **Facilitate PCR reaction set-up**
- **Reduced risk of contamination**
- **Direct gel loading**



## LongNova DNA Polymerase

**LongNova** DNA Polymerase is a mixture of *Pwo* and *Taq* polymerases. The thermostable **LongNova** DNA Polymerase catalyses DNA synthesis in a 5'→3' directions and shows 3'→5' exonuclease activity. It is characterized by high processivity and proofreading properties. It is ideal for long range PCR amplifications – up to 20 kb.

### Features:

- Wide range of product sizes – from 2 to 20 kb
- High proofreading properties (3'→5' exonuclease activity)

### Applications:

- Long PCR products
- Molecular cloning
- Site-directed mutagenesis and other methods which require high fidelity

Ordering information			
PRODUCT	CONCENTRATION	VOLUME	CAT. NO.
<b>LongNova</b> DNA Polymerase	2 U/μl	200 U	RP281
		1000 U	RP282
<b>LongNova-RED</b> DNA Polymerase	1 U/μl	200 U	RP281R
		1000 U	RP282R
<b>2x PCR</b> <b>LongNova-RED</b>	0.04 U/μl	100 reactions	RP85L
		1000 reactions	RP85L-10

## PCR Enhancers

### PCR Anti-inhibitor

The *PCR Anti-inhibitor* is a carefully composed mixture of alkaline proteins counteractive to various substances inhibiting the PCR reaction. Addition of the *PCR Anti-inhibitor* to a reaction mixture is an ideal way to eliminate the inhibitors derived from isolation process of a DNA used as a template for PCR. The *PCR Anti-inhibitor* should be added in 1:50 volume ratio to the PCR mixture.

Ordering information		
PRODUCT	VOLUME	CAT. NO.
<b>PCR Anti-inhibitor</b>	100 reactions	RP50
	500 reactions	RP51
<b>TaqSSB Protein</b>	50 μg	RP30
	250 μg	RP305
<b>5x GC-Additive</b>	1 ml	RP
	5x 1ml	RP



### Application

The *PCR Anti-inhibitor* is recommended for PCR with so called difficult DNA templates isolated from the specimen such as urine, saliva, sputum, blood, cell swabs, cerebrospinal fluid, biopsy specimen etc.

**M** – DNA Ladder M100-500      **K-** – PCR negative control

**S1 - S5** – amplification results obtained for reaction mixtures containing following SDS (PCR inhibitor) concentrations: S1: 0.1%; S2: 0.04%; S3: 0.02%; S4: 0.01%; S5: 0.005%

**K+** – PCR positive control (mixture does not contain SDS)

**A1 - A2** – amplification results obtained for mixtures containing 0.02% SDS and the following quantities of **PCR Anti-inhibitor** (reaction volume 50 μl): A1: 5 μl; A2: 1 μl

**A3 - A4** – amplification results obtained for mixtures containing 0.04% SDS and the following quantities of **PCR Anti-inhibitor** (reaction volume 50 μl): A3: 5 μl; A4: 1 μl

We own a substantial collection of thermostable single-stranded DNA-binding (SSB) proteins – PCR enhancers

## TaqSSB protein

A thermostable protein isolated from *Thermus aquaticus* binding single stranded DNA.

### Applications:

- Prevents PCR inhibition
- Increases amplification efficiency
- Increase selectivity and specificity of multiplex PCR
- Protects single stranded DNA from degradation
- Reduces secondary structures formation, which inhibit PCR
- Enhances amplification of difficult templates (e.g. rich in GC)
- Stimulates fidelity and processivity of *Taq* polymerase
- Reacts with RNA allowing size increase of synthesized cDNA in a RT-PCR reaction
- Stabilises ssDNA in the site-specific mutagenesis
- Facilitate in obtaining complete DNA digestion by restriction endonucleases

**Note!** We own a substantial collection of the SSB proteins. Remaining SSB proteins (from *Thermus thermophiles*, *Deinococcus radiopugnans*, *D. geothermalis*, *D. murrayi*, *Thermotoga neapolitana* and *T. maritima*) are produced upon request.

## dNTPs

Buffered solutions of ultrapure deoxyribonucleotides (dATP, dCTP, dGTP, dTTP) in various concentrations are ready to use in PCR, RT-PCR and Real-Time PCR reactions. The deoxyribonucleotides are a key component of an amplification reaction and their chemical purity is vital for the satisfactory results. They are supplied as lithium salts, which ensure prolonged stability and greater resistance to repeated freezing and thawing cycles. Thorough mixing of all 4 dNTPs guarantees high amplification efficiency and prevents incorrect incorporation of dNTPs in the synthesized strand.

### Properties:

- Highest purity appointed by an enzymatic synthesis and a high resolution HPLC (>99% dNTP, <0.9% dNDP)
- Increased sensitivity (allows template detection even in a very low copy number)
- High yield (DNase-, RNase-, ATPase- and pyrophosphatase-free)
- Effective and efficient PCR amplification even of very long amplicons (>10 kbp)
- Lithium salt dNTP guarantees high stability and greater resistance to repeated freezing and thawing cycles
- Free of dNTPs with modified bases and tetra- and pyrophosphates impurities (common PCR inhibitors)

## 5x GC-Additive

The PCR amplification of GC-rich DNA is often problematic due to the stable secondary structures of the DNA, which have high melting temperatures. These secondary structures interrupt DNA polymerase continuous movement along the DNA strand during the polymerization process, resulting in incomplete and nonspecific amplification. Many different methods and additives have been developed to facilitate template denaturation. The *5x GC-Additive* is an ideal tool for amplification on GC-rich templates. It changes the DNA behaviour upon heating and can be used with primer template – pairs with GC-rich content that do not work well with standard PCR conditions. The *5x GC-Additive* reduces the number of secondary structures and enables specific hybridization of primers. It is recommended to be used with *TaqNova* DNA Polymerase to obtain the most satisfactory results.

- **GC-rich templates**
- **High specificity**
- **Ideal for problematic templates, which fail with standard *Taq* DNA polymerases**

Ordering information			
PRODUCT	CONCENTRATION	VOLUME	CAT. NO.
<b>Equimolar mixtures of dATP, dCTP, dGTP, dTTP</b>			
dNTPs MIX 8 mM Total	2 mM of each dNTP	1 ml	RP61
dNTPs MIX 10 mM Total	2,5 mM of each dNTP	1 ml	RP63
dNTPs MIX 40 mM Total	10 mM of each dNTP	1 ml	RP64
dNTPs MIX 100 mM Total	25 mM of each dNTP	1 ml	RP65
<b>Sets of separate dATP, dCTP, dGTP and dTTP solutions</b>			
dNTPs SET 10 mM	10 mM of each dNTP (separately)	4x 0,2 ml	RP66
dNTPs SET 10 mM	10 mM of each dNTP (separately)	4x 1 ml	RP665
dNTPs SET 100 mM	100 mM of each dNTP	4x 0,2 ml	RP67
dNTPs SET 100 mM	100 mM of each dNTP	4x 1 ml	RP675



**Blirt S.A., DNA-Gdańsk**

Trzy Lipy 3/1.38, 80-172 Gdańsk, Poland  
[www.dnagdansk.com](http://www.dnagdansk.com) | [info@dnagdansk.com](mailto:info@dnagdansk.com)  
T: +48 58 739 61 50 | F: +48 58 739 61 51