



VIAFLO

Lightweight Electronic Pipettes

Your Pipette Tips Will Never Loosen, Leak Or Fall Off Again!

VIAFLO Lightweight electronic pipettes

VIAFLO electronic pipettes together with INTEGRA's wide range of GripTips form a perfect pipetting system. GripTips snap into place with minimal tip loading effort, providing a secure connection. GripTips never fall off and are always perfectly aligned, resulting in superior accuracy and precision results.

The ultra lightweight, well balanced design of the VIAFLO pipettes feature unsurpassed multichannel pipetting ergonomics, even during prolonged pipetting sessions.


The color screen together with the easy to navigate touch wheel user interface promote ease of set up and operation. All VIAFLO pipettes offer a choice of predefined pipetting modes covering typical applications as well as custom multi-step programs for specific tasks.



GripTips snap on

Do your pipette tips ever loosen, leak or fall off? This is a common issue in laboratories, caused by the use of universal pipette tips. Such tips require "hammering on" which stretches the shape of the pipette tip rim. The stretched pipette tip rim immediately starts moving back to its original shape, away from the tip cone. This can cause leaking, misaligned tips and even pipette tips that completely fall off the pipette!

The rim of the GripTips snap over the lobes, attaching firmly. As a result, GripTips will never accidentally fall off, even after several touch-offs. A shoulder provides a positive stop to prevent over-tightening of the tip. It can only be attached up to that point and no further. This means, the tip is either on or off, but cannot be in between. Unlike competitive pipettes, hammering the pipette into the tips to make sure they attach correctly is not required.



Versatile user interface

A color screen and proven touch wheel interface offer quick and comprehensive menu navigation. Simply by moving a finger over the touch wheel, menus and settings can be accessed. All menus are spelled out on the clear to read screen making memorizing abbreviations unnecessary.



A choice of ten predefined pipetting modes provides users with convenient solutions for the most common liquid handling tasks. Standard pipetting routines, such as repeat dispensing or serial dilutions can be executed with minimal set up time.

For more elaborate pipetting routines, up to 40 personalized custom programs can be created and saved. Custom programs can contain up to 98 individual steps based upon the following basic operations: Aspirate, Dispense, Mix, Purge and Prompt.


Ergonomic design

The new VIAFLO pipette offers perfect balance, an ultra lightweight design, unsurpassed operational comfort and is ideal for both left and right handed users. The effortless tip attachment and ejection forces add to the ergonomic benefits that accompany the use of VIAFLO electronic pipettes.



Repeat Dispense

Multiple dispensing of aliquots of the same volume without refilling the pipette tips after each dispense for rapid microplate filling.



Serial Dilute

Sequential transfer of a defined volume from one column to the next followed by a defined mix for convenient serial dilutions.



Manual Pipet

Only aspirate and dispense while the RUN button is being held down. Unlike a manual pipette, both volume and speed values can be set. Manual pipet mode can also be used to measure unknown volumes and perform titrations. The actual volume in the pipette tip is also shown on the display!



ASSIST / ASSIST PLUS

Automating multichannel pipettes for best results and unmatched ergonomics

The ASSIST pipetting robot family automates INTEGRA’s electronic multichannel pipettes to free scientists from routine pipetting tasks, such as serial dilutions, sample reformatting and reagent addition. It provides excellent pipetting reproducibility and sets a new standard for ergonomics.


ASSIST / ASSIST PLUS Ordering Information

	Part No.
ASSIST base unit	4500
ASSIST PLUS base unit	4505
Bluetooth module for VIAFLO / VOYAGER pipettes	4221



ASSIST PLUS

Ordering information

Instrument	Channels	Volume range	Part No.
	Single	0.5 - 12.5 µl	4011
		2 - 50 µl	4016
		5 - 125 µl	4012
		10 - 300 µl	4013
		50 - 1250 µl	4014
		100 - 5000 µl	4015
	Eight	0.5 - 12.5 µl	4621
		2 - 50 µl	4626
		5 - 125 µl	4622
		10 - 300 µl	4623
		50 - 1250 µl	4624
	Twelve	0.5 - 12.5 µl	4631
		2 - 50 µl	4636
		5 - 125 µl	4632
		10 - 300 µl	4633
	Sixteen	0.5 - 12.5 µl	4641
		2 - 50 µl	4646
		5 - 125 µl	4642

Instruments do not include a power supply. Choose your desired charging option below.

Charging options	Part No.
Linear stand (Holds up to 4 charging stations; mains adapter and charging stations sold separate)	3215
Mains adapter, for use with linear stand	3216
Charging station for linear stand (incl. connection cable)	3217
Charging/communication station for linear stand (incl. connection cable and USB cable)	3218
Mains adapter for electronic pipette	4200
Single pipette charging stand, incl. mains adapter	4210
Charging/communication stand, incl. mains adapter	4211
Carousel charging stand for 4 pipettes, incl. mains adapter	4215

GripTips for use with VIAFLO pipettes

For more GripTip options, including Low Retention, Wide Bore, and long/short versions, please refer to our website.

GripTips	Part No.	Part No.	Part No.
GripTip Size	Description	Sterile Filter	Sterile Non-Sterile
12.5 µl	5 Racks of 384 GripTips	4415	4414 4413
50 µl / 125 µl	5 Racks of 384 GripTips	4425	4424 4423
300 µl	5 Racks of 96 GripTips	4435	4434 4433
1250 µl	5 Racks of 96 GripTips	4445	4444 4443
5000 µl	5 Racks of 48 GripTips	4455	4454 4453

Bulk GripTips	Part No.	Part No.
GripTip Size	Description	Sterile Non-Sterile
12.5 µl	1 Bag of 1000 GripTips	4411
50 µl / 125 µl	1 Bag of 1000 GripTips	4421
300 µl	1 Bag of 1000 GripTips	4431
1250 µl	1 Bag of 500 GripTips	4441
5000 µl	1 Bag of 250 GripTips	4451
5000 µl	1 Bag of 100 GripTips, individually wrapped	4456

GripTip GREEN CHOICE Inserts	Part No.	Part No.
GripTip Size	Description	Pre-Sterilized Non-Sterile
12.5 µl	5 Inserts of 384 GripTips	4416 4412
50 µl / 125 µl	5 Inserts of 384 GripTips	4426 4422
300 µl	5 Inserts of 96 GripTips	4436 4432
1250 µl	5 Inserts of 96 GripTips	4446 4442

Charging options

Pipettes can be charged directly with the pipette mains adapter (4200). The single pipette charging stand (4210) and the carousel charging stand (4215) offer space saving charging options. The charging/communication stand (4211) additionally enables the communication of a pipette with a PC and VIALINK, a pipette management software from INTEGRA.



4210/4211

4215

The linear stand (3215) accommodates all types of INTEGRA pipettes. It can be easily converted into a charging stand, by equipping up to four charging stations (3217) and a mains adapter (3216).

Alternatively, a charging communication station (3218) can be equipped to enable the communication of a pipette with a PC and VIALINK.



3215

3217

3218

Reagent reservoirs

Disposable reservoir inserts designed to save you money and space, while keeping dead volume to a minimum. Available in 10 ml, 25 ml, and 100 ml configurations, as well as polystyrene or polypropylene and divided versions.



INTEGRA Biosciences AG
7205 Zizers, Switzerland
T +41 81 286 95 30
F +41 81 286 95 33
info@integra-biosciences.com

INTEGRA Biosciences Corp.
Hudson, NH 03051, USA
T +1 603 578 5800
F +1 603 577 5529
info-us@integra-biosciences.com

INTEGRA Biosciences Ltd.
Thatcham, Berks RG19 4EP, UK
T: +44 1635 797000
F: +44 1635 797001
info-uk@integra-biosciences.com



www.integra-biosciences.com

INTEGRA