


SERVAGel™ Neutral Gels


- Separation range from 5 kDa up to 200 kDa
- Long shelf life (> 1 year)
- Available as homogeneous and gradient gel
- Usage of different buffer systems



SERVAGel™ Neutral	SERVAGel™ Neutral	SERVAGel™ Neutral Gradient	SERVAGel™ Neutral Gradient
10 gels	10 gels	10 gels	10 gels
10 wells	12 wells	10 wells	12 wells
43222.01	43220.01	43223.01	43221.01

SERVAGel™ Neutral Gel HSE


- Optimized for High Speed Electrophoresis (HSE)
- 1D SDS PAGE: 20 minutes
- 2 D SDS PAGE: 30 minutes (2nd dimension)
- Available with 10 or 12 wells, 2D well



SERVAGel™ Neutral HSE	SERVAGel™ Neutral HSE	SERVAGel™ Neutral HSE 2D
10 gels	10 gels	10 gels
10 wells	12 wells	2D well
43246.01	43245.01	43247.01

SERVAGel™ Native Gels


- Suitable for both Blue and Clear Native PAGE
- Available as gradient gels 3–12 % and 4–16 %
- Available with 10 or 12 wells
- Ready-to-use anode, cathode and sample buffers are available



SERVAGel™ Native 3–12 %	SERVAGel™ Native 3–12 %	SERVAGel™ Native 4–16 %	SERVAGel™ Native 4–16 %
10 gels	10 gels	10 gels	10 gels
10 wells	12 wells	10 wells	12 wells
43251.01	43250.01	43252.01	43253.01

SERVAGel™ IEF Gels 2D Gels

- For vertical IEF in mini slab gel format
- Suitable for both native and denaturing conditions
- Standard IEF by cathodic sample application
- NEPHGE by anodic sample application for very basic proteins




SERVAGel™ IEF 3–10	SERVAGel™ IEF 3–10
10 gels	10 gels
10 wells	12 wells
43242.01	43240.01



SERVAGel™ Precast Gels
Your best choice in protein gel electrophoresis

- Precast polyacrylamide gels in mini format (10 x 10 cm)
- SDS free for native and SDS PAGE

SERVA
Electrophoresis



The SERVA PRiME™ System

PRiME™ - Premium Resolution in Minigel Electrophoresis

- BlueVertical™ PRiME™ mini gel electrophoresis apparatus
- SERVAGel™ PRiME™ precast mini gels at 10 x 10 x 0.7 cm
- High resolution in protein electrophoresis
- Separation range down to 3 kDa and up to over 200 kDa
- Run gels with Laemmli loading and running buffer
- Separation pattern comparable to standard Laemmli gels
- Long shelf life (> 1 year), ideal suited for stocking of gels in the lab
- Available as homogeneous and gradient gels
- With 10 or 12 sample wells, 2D well
- Running time approx. 35 minutes




It's PRiME™ Time

SERVAGel™ TG PRiME™

Homogeneous Gels

- Separation range from 3 kDa up to 200 kDa
- Run gels with Laemmli buffer system
- Available in 4 different AA concentrations: 8 %, 10 %, 12 % and 14 %
- Available with 10 or 12 sample wells




SERVAGel™ TG PRiME 8 %		SERVAGel™ TG PRiME 10 %		SERVAGel™ TG PRiME 12 %		SERVAGel™ TG PRiME 14 %	
10 gels	10 gels	10 gels	10 gels	10 gels	10 gels	10 gels	10 gels
10 wells	12 wells	10 wells	12 wells	10 wells	12 wells	10 wells	12 wells
43261.01	43260.01	43264.01	43263.01	43267.01	43266.01	43270.01	43269.01

SERVAGel™ TG PRiME™

Gradient Gels

- Separation range from 3 kDa up to 200 kDa
- Run gels with Laemmli buffer system
- Available in 3 different AA concentrations: 4–12 %, 4–20 % and 8–16 %
- Available with 10 or 12 sample wells




SERVAGel™ TG PRiME 4–12 %		SERVAGel™ TG PRiME 4–20 %		SERVAGel™ TG PRiME 8–16 %	
10 gels	10 gels	10 gels	10 gels	10 gels	10 gels
10 wells	12 wells	10 wells	12 wells	10 wells	12 wells
43274.01	43273.01	43277.01	43276.01	43280.01	43279.01

SERVAGel™ TG PRiME™

2D Gels

- Separation range from 3 kDa up to 200 kDa
- Run gels with Laemmli buffer system
- Available in 3 different AA concentrations: 12 %, 14 % and 8–16 %
- 2D well optimized for protein transfer from strip to gel



SERVAGel™ TG PRiME 12 % 2D	SERVAGel™ TG PRiME 14 % 2D	SERVAGel™ TG PRiME 8–16 % 2D
10 gels	10 gels	10 gels
2D well	2D well	2D well
43268.01	43271.01	43281.01

Equipment for PRiME™



BlueVertical™ PRiME™

- Mini electrophoresis tank system
- Ideally suited for SERVAGel™ PRiME™ gels
- Runs all gels with an outer format of 10 x 10 x 0.7 cm
- Unique, leak free clamp system
- Smart design - made in Germany

		BlueVertical™ PRiME
		1 Unit
		BV 104

BluePower™ 500x4 Power Supply

- 500 V, 1000 mA, 200 W
- Programmable (9 programmes, 9 steps each)
- For mini gel electrophoresis and blotting
- 4 x 2 outlets

		BluePower™ 500x4 Power Supply
		1 Unit
		BP-500x4

BluePower™ Supply Control Kit

- Monitoring of V, mA and W of each individual run
- Loading, storage and documentation of multistep settings
- Can be used with all SERVA BluePower™ supplies
- Compatible with all Windows® platforms (Windows® 98 and higher)

		BluePower™ Supply Control Kit
		1 Kit
		BP-PCSV01

SERVA
Electrophoresis

SERVA Electrophoresis GmbH
Carl-Benz-Str. 7
Fon: +49 6221 13840-0
info@serva.de
D-69115 Heidelberg
Fax: +49 6221 13840-10
www.serva.de

version 13/03