

Visit our Application Database to find
the right homogenizing solution for
your specific sample types!

www.omni-inc.com/apps

The Omni Bead Ruptor 24 is ideally suited for:

- DNA / RNA electrophoresis
- Liquid chromatography
- Northern, Southern, Western blot analysis
- Enzyme isolation
- Protein purification and production
- Food safety, GMO detection
- Cosmetic preparation
- Environmental surveys
- BSE Screenings
- Drug detections
- qPCR, RT-PCR
- And many more....



Specifications:

Processing Range: 25 μ L to 5mL

Power Rating: 100-230 VAC/ 50-60 Hz

Dimensions: 15.5" (W) x 17" (D) x 13" (H) [24" tall with lid open]

Weight: 24 kg (54 lbs.)

Compliance: CE Compliance

Run Time: 1 second—9:59 minutes

Number of Programs: 99 Storable Programs

Number of Cycles: 1—10

Performance Range: 0.8m/s—8m/s in increments of 0.15m/s

Dwell/Pause: 0—9:59 minutes



Biolabproducts GmbH
Hinter der Bahn 6
D-21376 G6denstorf
Tel.: 040 / 2000 4003
Fax: 040 / 2000 4001
email: info@biolabproducts.de
web: <http://www.biolabproducts.de>